# Fotoplast® Gel

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# Fotoplast® Gel

## Directions for use

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

### Product description:

Fotoplast® Gel is a light curing 1-component resin in a gel consistency suitable for modelling. It is specially suitable for repairs, for the modelling of vents and similar works at hard Fotoplast® or acrylic earmolds.

### Assembly of the syringe

### Processing:

For repairs apply Fotoplast® Gel onto a roughened and fat-free surface and cure in the polymerization light unit. In order to apply the venting clean and decrease the inside of the shell. Place the venting wire and surround it by Fotoplast® Gel and cure it in a suitable UVA light unit.

**Polymerization:** In the Polyflux with 300-400 nm  
**Curing time:** approx. 2 min.  
**Storage:** 18 - 28 °C, keep tube tightly closed  
**Available in:**

10 ml clear transparent	REF 44811	20 ml clear transparent	REF 44691
10 ml opaque H (yellowish)	REF 44821	20 ml red transparent	REF 44612
10 ml red transparent	REF 44831	20 ml blue transparent	REF 44622
10 ml blue transparent	REF 44841	80 ml clear transparent	REF 44791
20 ml opaque H (yellowish)	REF 44651		

### Fotoplast® Gel transparent:

**Hazard statements**  
**H360Df** May damage the unborn child. Suspected of damaging fertility.

### Precautionary statements

**P280** Wear protective gloves/protective clothing/eye protection/face protection. / **P201** Obtain special instructions before use. / **P202** Do not handle until all safety precautions have been read and understood. / **P308+P313** IF exposed or concerned: Get medical advice/attention. / **P405** Store locked up. / **P501** Dispose of contents/container in accordance with local/regional/national/ international regulations.

### Additional information:

Contains Diphenyl(2,4,6-trimethylbenzoyl)phosphinoxide. May contain an allergic reaction.

**Contains:** tetrahydrofurfuryl-2-methacrylate

### Fotoplast® Gel opaque:

**Hazard statements**  
**H315** Causes skin irritation. / **H319** Causes serious eye irritation. / **H317** May cause an allergic skin reaction. / **H335** May cause respiratory irritation.

### Precautionary statements

**P261** Avoid breathing dust/fume/gas/mist/vapours/spray. / **P280** Wear protective gloves/protective clothing/eye protection/face protection. / **P305+P351+P338** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. / **P304+P340** IF INHALED: Remove person to fresh air and keep comfortable for breathing. / **P405** Store locked up. / **P501** Dispose of contents/container in accordance with local/regional/national/international regulations.

### Other hazards

Uncontrollable exothermic product polymerisation possible in case of storage above 60 °C or exposure to direct sunlight.

**Contains:** hydroxypropylmethacrylate, monoester with propane-1,2-diol, Acrylic resin, 2-hydroxypropyl methacrylate, 2-hydroxyethyl methacrylate

Dreve Otoplastik GmbH is not liable for any damages caused by improper application of the material.

2017-12-01 / Rev. 1  
 Dreve Otoplastik GmbH  
 Max-Planck-Str. 31 • 59423 Unna / Germany  
 Tel.: + 49 2303 8807-0 • Fax: + 49 2303 82909  
 E-Mail: otoplastik@dreve.de • Internet: www.dreve.com