
Screening Audiometer

QUICK USER'S HANDBOOK

resonance

CE
0051

resonance

Release history

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this page is intentionally left blank

Table of contents

1-Introduction and intended use	29
2-Description	30
2.1 Precautions	30
2.2 Warranty and maintenance	31
2.3 Unpacking and checking the R07A	32
2.4 Accessories	32
2.5 Connecting the accessories	34
2.6 Batteries replacement	35
2.7 Symbols present on device	36
3-Function of buttons and display	37
3.1 Power On and Off the device	37
3.2 Pure Tone presentation	38
3.3 Function buttons	38
3.4 Display	39
4-Performing tests	40
4.1 Pure Tone Audiometry (Air Conduction)	40
4.2 Hugson-Westlake Auto-threshold test (Air Conduction)	41
5-Menu	42
5.1 Patient	42
5.2 Setup menu	43
5.3 Frequencies	44
5.4 Calibration	44
5.5 About	44
6-Care and Maintenance	45
7-Safety Regulations	45
7.1 Electrical safety	45
7.2 Device handling and operation	46
7.3 Measurement safety	46
8-Technical data	46
9-Electromagnetic compatibility	48

1. Introduction and intended use

RESONANCE thanks you for having purchased one of its medical devices.

This manual contains information and warnings aiming at ensuring the correct and safe use of the R07A Screening audiometer part number 4300104200.

Carefully store this manual away in a safe place for each future consultation.

The R07A Screening audiometer is designed and manufactured to meet or exceed all quality and safety requirements, and has been certified with CE-symbol according to Medical Directive. All the functions are software-controlled, upgrading the software at a later date will be simple and cost-effective.

The R07AScreening audiometer is intended to be used by an audiologist, hearing healthcare professionals, or trained technicians in a quiet environment.

The user must therefore have an adequate knowledge of audiological diagnostic equipment, expertise on the examination methods used and knowledge of English.

The instrument is intended for use in quiet environments, ideally sound treated to <35 dB (A).

It is recommended that the device be operated within an ambient temperature range of 15 - 35 degree Celsius.

The R07A Screening audiometer is designed to be a device for screening for hearing loss; in order to assess reliable audiometric results, some variables factors must be kept in mind as they may influence and/or jeopardize test results as:

Test room installation

- noise level in the test room
- correct calibration of the audiometer

Audiometric tests

- correct position of the transducers
- full instructions to the patient to obtain reliable responses
- testing methods, ascending or descending

Failing to account for these factors may lead to inaccurate results.

2. Description

2.1 Precautions

WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.

NOTICE is used to address practices not related to personal injury.

- To maintain the safety and full functionality of the R07A, it is important that instructions set forth in this manual are strictly followed.
- Ensure compliance with environmental conditions (during transport, storage and operating). If the device is moved from a cold location to a warmer one or vice versa, there will be a risk of condensation build up. If condensation occurs the device must be allowed to achieve normal temperature before it is turned on.
- To avoid damage to the device and leakage from the batteries, batteries should be removed whilst the device is being stored.
- R07A does not have specific ingress protection. Care should be taken to avoid liquid ingress. Immersion of the device may result in serious damage.
- Although the instrument fulfils the relevant EMC requirements, precautions should be taken to avoid unnecessary exposure to strong electro-magnetic fields, e.g. from mobile phones, microwaves, radiotherapy devices, etc.; this may induce malfunctions in their performance. If the device is used adjacent to other equipment it must be observed that no mutual disturbance appears.
- The R07A is intended to be powered from batteries or from a medical CE power supply offered as an option. Using other supplies is not allowed. When used, position it so that the power plug can be easily disconnected.
- The device should be placed at least 20 cm from a heat source.
- Be sure to use only stimulation levels, which will be acceptable for the patient.
- The transducers supplied with the R07A are calibrated to this device (see DHR documents supplied with the device); exchange of transducers will require a recalibration.
- Sudden impacts or drops may cause damage to the device, which may affect its performance and safety. Repairs should only be undertaken by a qualified technician.

If the R07A is connected to other devices with their own power source, the system conforms to EN 60601-1 must be checked under the responsibility of the installer; in the event that isolation devices are used to achieve such conformity, these must comply with EN 60601-1. Any person who connects external equipment to the device connectors has formed a system and is therefore responsible for the system to comply with the requirements of IEC 60601-1.

If in doubt, contact your RESONANCE technician or local distributor for help.

Device not to be used in areas in which flammable anaesthetics or flammable agents for disinfection or skin cleaning are applied.

Electric and electronic waste may contain hazardous substances and therefore has to be disposed of separately.

Such products will be marked with the crossed-out wheeled bin shown.

The cooperation of the user is important in order to ensure a high level of reuse and recycling of electric and electronic waste.

Failure to recycle such waste products in an appropriate way may endanger the environment and consequently the health of human beings.

Disposal of batteries must be made in accordance to national regulations.

2.2 Warranty and maintenance

RESONANCE warrants that:

- The R07A is free from defects in material and workmanship under normal use and service for a period of 36 months from the date of device's installation by authorized personnel to the first purchaser.
- Accessories are free from defects in material and workmanship under normal use and service for a period of 12 months from the date of device's installation by authorized personnel to the first purchaser.

Furthermore, this warranty shall not apply to, and M.R.S. shall not be responsible for, any loss arising in connection with the purchase or use of any RESONANCE device that has been:

- repaired by anyone other than an authorized RESONANCE service representative;
- altered in any way so as, in RESONANCE judgment, to affect its stability or reliability;
- subject to misuse or negligence or accident, or which has had the serial number altered, effaced or removed; or improperly maintained; or used in any manner other than in accordance with the instructions furnished by RESONANCE.
- This warranty is in lieu of all other warranties, express or implied, and of all other obligations or liabilities of M.R.S. , and M.R.S. does not give or grant, directly or indirectly, the authority to any representative or other person to assume on behalf of M.R.S. any other liability in connection with the sale of RESONANCE products.
- The R07A may be repaired only by your dealer or by a service company authorised by RESONANCE, recommended by your dealer.
- Only original RESONANCE spare, accessories, cables or consumable parts must be used for service or repair the R07A.

In order to ensure that your device works correctly, the R07A audiometer should be checked and calibrated at least once a year. We recommend to send the R07A to your RESONANCE dealer for this service, including the headphone and all of the accessories. Routine maintenance includes checking the status of the transducer cables, of medical power supply cables (if present), headset cushions, verifying the integrity of the cabinet of the instrument. The verification of instrument calibration must be conducted by qualified service personnel authorized by the manufacturer via appropriate instrumentation. This must provide the level control and the frequency of all signals emitted by the transducers supplied with the instrument.

In case of replacement of transducers and / or parts, use only parts provided by the manufacturer. In the case of replacement of transducers, a new calibration of the instrument is required.

2.3 Unpacking and checking the R07A

On checking for packaging and content damage immediately after the device's reception; should any sign of damage or tampering be noted, please notify the carrier immediately.

If the device has been damaged during transportation, the device should be checked for any electrical or mechanical defects.

Please contact your RESONANCE dealer in case of any damage is identified.

Keep all original packaging for future use.

The identification label on the rear plate holds the serial number. This should be checked and written down for later service claims.

2.4 Accessories

List of accessories

HB-7 Audiometric headset with DD45 transducers

Patient response pushbutton

100 sheets audiogram charts

Blue/Red pencil

Multilingual CE Operation manual

MDS Software

USB pen drive including: DHR documents, MDS software installer, Multilingual CE Operation manual

Carrying bag

4 AA 1,5 V batteries

Optionals available

HB-7 Audiometric headset with TDH39 transducers

IP30 Insert phones

Medical CE approved power supply

ADC Audiocups Noise Excluders

Available tests

Pure Tone,

Automatic audiometry based on Hughson-Westlake algorithm

Image	Description	Standard / Optional	Plug
	DD45 headphone*	●	ACR/ACL (A-B)
	IP30 insert headphone*	○	ACR/ACL (A-B)
	Patient pushbutton*	●	PAT RESP (C)
	Medical power supply: Model: MENB1010A0503(xyz) (xyz) = B01 changeable direct plug-in (xyz) = Q01 with power cord (xyz) = N01 with power cord	○	5V DC (D)
	As alternative: Medical CE approved power supply: Model: ME10A0503(xyz) (xyz) = B01 changeable direct plug-in (xyz) = Q01 with power cord (xyz) = N01 with power cord	●	5 VDC (1)

● = standard

○ = optional

* = applied part

The applied parts are represented by the transducers ad pushbutton which are placed in contact with the patient. They are indicated by an asterisk (*).

2.5 Connecting the accessories

All the connections jacks can be found on the left side of the R07A.

Please connect all the cables before powering on the device.


Figure 1 – side view of the R07A

Reference	Symbol	Function
A	ACR (Red)	Socket for AC Right headphone channel
B	ACL (Blue)	Socket for AC Left headphone channel
C	Pat. Resp. (Black)	Socket for patient pushbutton
D	5V DC	Socket for external power supply
E		Mini USB port for PC communication
F		Battery compartment


Figure 2 – rear view of the R07A

2.6 Batteries replacement

Remove the cover of the battery compartment located on the bottom shell.

Remove the old batteries and insert the new ones, matching the polarity as detailed in the label present of the compartment.

To avoid damages, insert each battery starting from the positive polarity.

See images here below for more details.

Positive polarity first


Respect polarity as in label

2.7 Symbols present on device

	Refer to instruction
	Refer to manual for transducer and power supply connection
	Class II device. (Symbol refers to external medical power supply with double insulation)
	BF type applied parts
	Electric and electronic waste may contain hazardous substances and therefore has to be collected separately
	Conformity marks to RoHS directive
	Together with name and address indicates the manufacturer of the instrument
	Conformity marks to directive 93/42/CEE – notification body 0051 (IMQ)
	Year of production of the device
	Slave USB plug

3. Function of buttons and display

3.1 Power on and off the device

To turn on the device, press the green button Main switch button.

To switch off the device, press the green button "Main switch" and hold for a few seconds until the message "POWER OFF!" is displayed on the screen.

Before powering on the device, ensure that all the accessories are already plugged-in; do not remove them during test execution, power off the device before to disconnect the cables.

The R07A can also power off automatically after 1, 2 or 5 minutes (or never) depending on the settings.

3.2 Pure Tone presentation

To send the tone stimulus to the patient, you must use these three commands:

Signal: Touch the silent Signal switch; because it is a soft-touch switch, there is no need to press on it.

Level: Rotate the left knob (dB) to increase or decrease the level of the stimulus.

Frequency: Rotate the right knob (Hz) to select different frequencies. Please note that only activated frequencies can be used. See point 5.2 on this manual.

3.3 Function buttons

R/L: With this button you can select which ear you deliver the stimulus. The display shows the selected setting.

INPUT: This button is used to select between steady, pulsed or warble tone presentation. You can select different modulation rates for the Pulsed tone; see point 5.2 on this manual.

On the bottom right side of the display a symbol will be shown to identify the selected tone presentation modality:

TALK FORWARD: Pressing this button will activate the talk forward communication with the patient. This option is useful when you have to instruct the patient. The display shows the talk forward level that can be adjusted rotating the one of the two knobs (Frequency or Level).

To exit from this function, press the Talk Forward button again.

ON-OFF: This is a multifunction switch. When the device is off, to turn on the audiometer, press on this switch. The same operation can be made when the device is on; to switch off the device keep pressed this switch for a few seconds.

To avoid loss of any data, be sure to store the latest tests results before switch off the device.

Once on the device, this switch allows to access main menu.

F1 - F2: These two soft keys are used to perform different actions. Depending from the active menu, the two soft keys can be used for exit, select or store the selected function.

3.4 Display

SIGNAL: A stimulus presentation indicator is provided in the top right corner on the display header.

RESPONSE: When using the patient response switch, a response is indicated in the middle of the display header.

BATTERY STATUS: The power status of the R07A is indicated in the top right corner of the display header.

Empty:

Full:

DC powered:


The icon will change depending on whether the device is powered via its external power supply or batteries.

When powered by batteries, the battery icon will change depending on the battery power level. When batteries are running low, the display will show a low battery and flash.

The automatic power off setting of the device can be adjusted to different time intervals; or set to never power off – please see Setup section for details.

Only medical power supply approved and supplied by the manufacturer should be used to power the device.

4. Performing tests

4.1 Pure Tone Audiometry (AC)

Pure Tone Audiometry represents the basic test of auditory investigation aiming at evaluating the auditory function and assessing the hearing threshold. The term hearing threshold indicates the minimum level perceived by the test subject.

The test is performed by administering pure tones using the included earphones (Air Conduction – AC). Threshold estimation can be done either for one ear and then for the other, for all tones, or on both sides for each tone.

The R07A can specify the AC loss but cannot distinguish between abnormality in the conductive mechanism and sensorineural mechanism.

The patient should sit at a distance of at least 1 m from the device.

Before perform the test, eliminate any obstructions which will interfere with the placement of the earphone cushions on the ear (e.g. hair, eyeglasses, rings).

Place the right headphone (identified by red color on the transducer) on the right side and the left (identified by blue color on the transducer) on the left.

Ensure that the cushions are clean to avoid any potential infection's transmission.

Adjust the headband of the headphones sitting on the two metallic bars, to achieve the correct pressure on the ear and the correct alignment of transducer to ear canal.

Instruct the patient, before performing the test, to signal pressing the pushbutton (or to raise his/her hand) when he/she hears a sound, and to keep the pushbutton depressed (or his/her hand up) as long as he/she hears it.

Threshold determination: The test begins at 1000 Hz on the patient's better ear.

Present pure tones at 1000, 2000, 4000, 8000 Hz then 500, 250 and 125 Hz*.

Start presenting the tone at a level loud enough to ensure audibility. The patient will thus be able to identify the type of signal to be expected exactly and give clear responses.

Should levels close to threshold be used when starting the test, the patient may fail to respond as he/she is unsure of what he/she should hear. Should the level prove insufficient to elicit clear responses, increase the level until the patient responds reliably.

To assess threshold level, proceed as follows:

- When the stimulus is clearly perceived by the patient, reduce progressively its level in 10 dB steps.
- When reaching a level, at which the patient fails to respond, further decrease the level by 10 dB and, interrupting the tone manually, begin increasing its level in 5 dB steps, until the patient responds again.
- Decrease this level once more until the patient fails to respond.

The true pure tone test threshold corresponds to the level evoking patient's response in 50% of presentations.

In tinnitus cases, use of pulsed tones (rather than continuous) is preferable.

Avoid any noisy environments because background noise can produce false test results, especially at lower frequencies.

* 125 Hz is not standard in the British Recommended Procedure

4.2 Hugson-Westlake Auto-threshold (AC)

The test, controlled by the patient's pushbutton according to the Hugson-Westlake method, provides a modified and abbreviated version of the ISO 8253-1 standard for recording automatic AC threshold without masking.

This test is especially used for industrial audiology and hearing conservation programs.

Hugson Westlake is an automatic pure tone test procedure. Threshold is defined as 2 out of 3 correct responses at a certain level in a 5 dB increase and 10 dB decrease test procedure. The Hugson Westlake method is used to obtain pure tone thresholds automatically.

How perform the test: Test patient must be instructed to keep his/her pushbutton pressed as long as he perceives the tone and release it as soon as the tone is no longer heard.

To activate the Auto-threshold test please refer to paragraph 5.2 Setup menu.

Once activated, entering the *Menu* you can select the test pressing F2.

The screen will display the test interface.

Familiarisation

Before starting the test, you can perform a familiarisation to train the patient.

Use the buttons as described on the paragraph 3.2 to act the familiarization process: present a tone at 1000 Hz which can easily be perceived (i.e. 50 dB) If necessary, increase with steps of 10 dB until the tone is clearly perceived.

If the familiarisation process is successfully completed, you can start the test pressing F1.

Threshold Determination

The hearing threshold is defined as the lowest level at which at least 50% of the stimuli elicit a response. This threshold is found by the following procedure.


- 1) Present a tone which is 10 dB lower than the level at which familiarisation was finished.
- 2) Decrease the level in steps of 10 dB until response fails.
- 3) Increase the level in steps of 5 dB until the subject responds again.
- 4) Repeat 2) and 3) two or three times until the threshold appears at the same level to at least 50% of presentations.

The time intervals in between the stimuli should be varied to prevent the patient predicting a pattern to the stimulations.

- 5) Change to the next frequency and repeat the procedure until all frequencies are measured. Repeat the procedure at 1000 Hz. If the difference to the previously found threshold is then 5dB or less proceed to test the other ear. If the difference is 10 dB or greater, repeat the test at the other frequencies, until agreement within 5 dB has been obtained.
- 6) Proceed until both ears have been tested.

Once completed the test, store the results pressing F2.

5. Menu

To access the R07A *Menu* press ON/OFF switch when you are on the Test execution modality. To browse up or down in the *Menu*, rotate either knob.

Entering this menu you can set or select:

Patient
Setup
Frequencies
Calibration (for trained personnel only)
About

5.1 Patient

To access the R07 Patient management, enter on *Menu* then select the associated function by pressing F2 soft key or pressing one of the two knobs.

Rotating the knob, you can toggle through a list of 50 patients available.

The patients are identified only by a progressive number.

Once a patient is selected, three options are available:

F1- DELETE The data stored to this patient will be deleted and removed from the internal database.
Please confirm if you want to delete data pressing F2 before exit.

F2- SELECT If the patient is in "FREE" condition, pressing on the F2 soft key, the test result record will be activated related to the selected patient. The patient ID number will be display on the test screen.

F2- VIEW If the patient is in "USED" condition, pressing the F2 soft key, all the stored results can be viewed. Rotating the knob all the frequencies can be viewed.

Pressing EXIT the device will exit to the previous menu without saving any new selection.

5.2 Setup menu

To access the R07A setup menu press F2 soft key or one of the two knobs.
To browse up or down in the Setup Menu, rotate either knob.

HEADSET

Press CHANGE to select an alternative transducer; DD45, Insert earphones or TDH39. The default setting is DD45. Press STORE to save the new setting or EXIT. On exiting, the setting will not be updated.

START EAR

Press CHANGE to toggle between RIGHT and LEFT ear as the default ear to start to perform a test. Press STORE to save the new setting or EXIT. On exiting, the setting will not be updated.

SIGNAL (Tone presentation modality)

Press CHANGE to toggle between NORMAL or REVERSE. NORMAL: Tone is presented as long as the signal switch is pressed. REVERSE: Tone will be interrupted if signal switch is activated. REVERSE mode results in higher battery power consumption. Press STORE to save the new setting or EXIT. On exiting, the setting will not be updated.

POWER OFF settings

Press CHANGE to toggle between 3, 4 or 5 minutes. If powered by batteries, the device will shut down after "Power-off time" as entered in the setting. The countdown of this time will start when the display backlight will be switched off taking count of the "LCD-off time" setting. Press STORE to save the new setting or EXIT. On exiting, the setting will not be updated.

With power supply connected, the device will NOT turn-off automatically.

STEP (Stimulus Level Steps)

Press CHANGE to select between 1, 2 or 5 dB steps. Press STORE to save the new setting or EXIT. On exiting, the setting will not be updated.

DEFAULT dB (Default Level)

Press CHANGE to choose between a range from 20 up to 60 dB in steps of 5 dB the default level when changing ear side. The default setting is 40 dB. Press STORE to save the new setting or EXIT. On exiting, the setting will not be updated.

PULSED (Pulse Length)

Press CHANGE to toggle between 0.5, 1 or 2 Hz to adjust the length of the pulsed tone. Press STORE to save the new setting or EXIT. On exiting, the setting will not be updated.

LCD OFF settings

Press CHANGE to toggle between 0, 30, 60, 90 or 120 seconds to adjust the activation time of "save batteries modality". In this condition, the display backlight is switched off to save the battery power. All the functions of the R07A are still active, acting on any button or soft key the display backlight will automatically resume. These settings are active only when the device is powered with batteries; if plugged to external power supply LCD OFF is disabled and backlight is always active.

CONTRAST

Press CHANGE to toggle between setting ranging from 0 (dark) up to 5 (bright) to adjust the LCD contrast. Press STORE to save the new setting or EXIT. On exiting, the setting will not be updated.

LANGUAGE

Press CHANGE to toggle between the languages available.

AUTOTHRS

Press CHANGE to toggle between NO (Auto-threshold test not active), QUICK or NORMAL. If active, on the top of the main menu will be displayed the possibility to enter on the Auto-threshold test mode. Press STORE to save the new setting or EXIT. On exiting, the setting will not be updated.

AUTOFREQ

Press CHANGE to enable or disable this function. If active, this function will switch automatically to the next frequency, based on the related setting, once a value is stored. This function will help the operator to reduce the testing time.

FACT RESET

This function resets the device to default settings.

Press CHANGE to select YES or NO, press STORE to confirm your selection.

5.3 Frequencies

To access the R07A Frequencies menu press F2 or one of the two knobs.

10 frequencies are available to enable or disable: 125, 250, 500, 750, 1500, 2000, 3000, 4000, 6000 and 8000 Hz.

1000 Hz cannot be disabled.

Once a frequency has been selected by rotating a knob, simply press on "OFF" or "ON" soft key to enable or disable the selected frequency.

Pressing the soft key "Exit", the new settings will be stored and the screen will show again the General menu.

5.4 Calibration

This section is for trained personnel only and allow to perform the annual calibration of the device.

To enter on this section you need a password.

5.5 About

Press F2 soft key or one of the two knobs to access the information in the About section.

6. Care and maintenance

Switch off the device and disconnect the power plug before cleaning!

Always disconnect the mains conductor during the cleaning process, and be careful that no fluid enters inside of the device or the accessories; no alcohol or spirits should be used. The instrument requires cleaning the housing with a damp cloth (do not use aggressive cleaning agents); absolutely avoid the penetration of liquids.

It is recommended that parts which are in direct contact with the patient (e.g. earphone cushions) are subjected to standard disinfecting procedure between patients.

After each test, it is necessary to disinfect the headset cushion before use on a new patient. This includes physically cleaning and use of a recognised disinfectant. Individual manufacturer's instruction should be followed for use of this disinfecting agent to provide an appropriate level of cleanliness.

In case of use R07A with insert earphones, never insert or in any way use the insert headset without a new clean non defect test tip supplied by RESONANCE.

If the surface of the device or parts of it are contaminated, it can be cleaned using a soft cloth moistened with a mild solution of water and dish washing cleaner or similar. The use of organic solvents and aromatic oils must be avoided all the times.

Before performing test:

- It is necessary to disinfect the headset cushion before use on a new patient.
- The tips of the Insertion tips are disposable: use them only once and replace them before using them with a new patient.

Daily verify:

- Presence of signals in the transducers provided
- Correct operation of all other controls provided

Weekly verify:

- Good condition of all accessories cables and related jack-plugs

Annually:

- Perform calibration and routine maintenance check

7. Safety regulations

7.1 Electrical Safety

The R07A audiometer is in compliance with Class BF of EN 60601-1, in any case, the device is not intended to be used in environments dealing with explosive material or equipment.

7.2 Device handling and operation

The device should be checked once daily, as per operations described on paragraph 6.

The device is designed to be handled and operated by trained personnel like Audiologist, ENT doctors or personnel with similar qualifications ONLY.

7.3 Measurement Safety

As a manufacturer, we recommend to perform an annual inspection and calibration in order to ensure safety and quality of the measurement.

The annual check-up can be performed only by RESONANCE distributor or authorised service center. M.R.S. will not be liable for any failure to comply with the specified inspection date, according to the Medical product law. The use of uncalibrated and uninspected audiometer is strictly prohibited.

8. Technical data

Reference standard:	EN 60645-1/ANSI S3.6, Type 4
Safety:	EN 60601-1
EMC:	EN 60601-1-2
Medical CE Marks:	The Manufacturer's Quality Management System complies to the Annex II of Medical Device Directive 93/42/EEC. The notification body IMQ, Identification No. 0051, has approved the Quality System.
Inputs:	Medical approved Power supply; Patient pushbutton
Outputs:	AC Left and Right
Attenuator:	-10 to 100 dB HL in 1, 2 or 5 dB steps
Tone presentation:	Manual or Reverse (selected in Setup menu); Pulsed at 0.5, 1, 2 Hz (selected in Setup menu); Manual On/Off; Pure tone or Warble
Talk Forward:	Built-in talk forward microphone 50 - 100 dB SPL adjustable
Auto Threshold:	Patient controlled Hughson Westlake procedure according to ISO 8253-1
Store function:	Soft key store button and internal memory for 50 patients and related tests results. Stored measurements can be viewed on the build in display.
Modulation:	Warble +/- 5% 5 Hz.

Available frequencies and maximum Intensities:

Maximum Intensities											
Frequency	125 Hz	250 Hz	500 Hz	750 Hz	1 KHz	1.5 KHz	2 KHz	3 KHz	4 KHz	6 KHz	8 KHz
DD45	70	90	100	100	100	100	100	100	100	100	90
IP30	80	90	100	100	100	100	100	100	100	95	80

Calibration

Air Conduction: ISO 389-1/ANSI S3.6 (DD45)

AC-Headphone DD45

Reference equivalent threshold according to ISO 389-1 with coupler IEC 60318-3 Force 4-5 N

Frequency (Hz)	dB
125	47.5
250	27.0
500	13.0
750	6.5
1000	6.0
1500	8.0
2000	8.0
3000	8.0
4000	9.0
6000	20.5
8000	12.0

AC-Insert Phones IP30

Reference equivalent threshold according to ISO 389-2 with coupler IEC 60318-5

Frequency (Hz)	dB
125	26
250	14
500	5.5
750	2
1000	0
1500	2
2000	3
3000	3.5
4000	5.5
6000	2
8000	0

General Dimensions

W x D x H: 240 x 190 x 47 mm

Weight: 1.0 Kg – including batteries and headset

Environmental conditions

Transport and storage: Temperature: from -20 to +50 °C

Humidity: from 10 to 95% without condensation

Operation: Temperature: from 10 to +35 °C

Humidity: from 30 to 90% without condensation

Power

Batteries: 4 type AA batteries / External Power supply

Automatic battery on/off switching

Battery status indication

External power supply: 5VDC 100 - 240 V, 50 - 60 Hz

Construction

Plastic ABS cabinet

9. Electromagnetic compatibility

Guidance and manufacturer's declaration - Electromagnetic EMISSION		
The R07A audiometer is intended for use in the electromagnetic environment specified below. The customer or the user of the R07A should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF Emissions CISPR11	Group 1	The R07A uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR11	Class B	The R07A is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies building used for domestic purposes.
Harmonic Emissions IEC/EN 61000-3-2	Complies	
Voltage fluctuations/flicker emission IEC 61000-3-3	Complies	

Producer's guide and declaration - Electromagnetic IMMUNITY			
The R07A is suitable for use in the specified electromagnetic environment. The purchaser or user of the R07A should assure that it is used in an electromagnetic environment as described below:			
Immunity test	IEC 60601-1-2 Test level	Compliance Level	Electromagnetic Environment
Electrostatic discharge (ESD) IEC 61000-4-2	8 kV contact 2/4/8/15 kV air	IEC 60601-1-2 Test level	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.

Radiated electromagnetic field IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	IEC 60601-1-2 Test level	Portable and mobile RF communications equipment should be used no closer to any part of the R07A including cables. Recommended instance: $d = 1,2\sqrt{P}$ 80 MHz to 800 MHz $d = 2,3\sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum power output of the transmitter in watt (W) based on the transmitter's producer and d is the recommended distance in meters (m).
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines > 3 m	IEC 60601-1-2 Test level	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	0.5/1 kV differential mode 0.5/1/2 kV common mode	IEC 60601-1-2 Test level	Mains power quality should be that of a typical commercial or hospital environment.
Conducted disturbances induced by RF fields IEC 61000-4-6	3 V 150 kHz to 80 MHz	IEC 60601-1-2 Test level	Portable and mobile RF communications equipment should be used no closer to any part of the R07A including cables. Recommended distance d = $1,2\sqrt{P}$
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% Un for 0.5 cycle 0 % Un for 1 cycle 70 % Un for 25 cycles 0 % Un for 5 s	IEC 60601-1-2 Test level	Mains power quality should be that of a typical commercial or hospital environment. If the user of the R07A requires continued operation during power mains interruptions. it is recommended that the R07A be powered from an uninterrupted power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	IEC 60601-1-2 Test level	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: At 80 MHz and 800 MHz, the highest frequency range is applied

Note: These guidelines may not be applicable in all situations. Electromagnetic propagation is influenced by the absorbance and deflection of structures, objects and people.

- a. The intensity of the fixed transmitter's field such as radiophone (cellular/cordless) bases and land mobiles, amateur radios, FM and AM radio and television, theoretically cannot be determined precisely. To evaluate the electromagnetic environment due to fixed RF transmitters, a local electromagnetic investigation should be considered. If the measured field intensity in which the R07A will be used, is above the applicable conformity level of RF described above, the R07A should be under observation to verify normal function. In case of anomalous function additional attention may be required, such as changing position or orientation of the R07A.
- b. Above the frequency range from 150 kHz to 80 MHz, the intensity of the field must be below 3 V / m.

Recommended distance between communication portable and mobile and R07A

R07A is for use in an electromagnetic environment in which irradiated RF interferences are controlled. The clients or users can minimize electromagnetic interference by maintaining a minimal distance between portable and mobile RF transmitters and R07A as indicated below, based on the max output power of the communication devices.

Max nominal power output of the transmitter W	Distance based on the frequency of transmitter m		
	150 KHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.2
100	12	12	23

For those transmitters with a max nominal output power not mentioned above , the advised distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the max output power in watt (W) according to the transmitter's builder.

Note 1: at 80 MHz and 800 MHz, the distance applied refers to the highest range frequency

Note 2: These guidelines may not be applicable in all situations. Electromagnetic propagation is influenced by the absorbance and deflection of structures, objects and people.

this page is intentionally left blank

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Image	Description	Standard / Optional	Plug
A pair of black over-ear headphones with a coiled cable.	DD45 Kopfhörer*	●	ACR/ACL (A-B)
A blue and red earpiece with a small microphone attached.	IP30 Einstekkopfhörer*	○	ACR/ACL (A-B)
A blue rectangular button with a small antenna-like protrusion.	Patientenantworttaste*	●	PAT RESP (C)
A collection of medical power supply components including a main unit, multiple international power adapters, and various cables.	Medizinische Stromversorgung: Model: MENB1010A0503(xyz) (xyz) = B01 austauschbarer direkt Kontaktstecker (xyz) = Q01 mit Anschlusskabel (xyz) = N01 mit Anschlusskabel	○	5 VDC (D)
A smaller black power supply unit with a short cable.	As alternative: Medical CE approved power supply: Model: ME10A0503(xyz) (xyz) = B01 changeable direct plug-in (xyz) = Q01 with power cord (xyz) = N01 with power cord	●	5 VDC (1)

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Imagen	Descripcion	Standard / Opcional	Plug
A pair of black over-the-ear headphones with a coiled cable.	Fonos DD45*	●	ACR/ACL (A-B)
A pair of blue earphones with a red and blue cable, designed for insertion into the ear.	Fonos de Inserción IP30*	○	ACR/ACL (A-B)
A black patient response button with a flexible cable.	Boton respuesta paciente*	●	PAT RESP (C)
A collection of power supply components including a main unit, multiple international power cords, and a changeable direct plug-in cord.	Fuente de Poder: Modelo: MENB1010A0503(xyz) (xyz) = B01 plug irectointercambiable (xyz) = Q01 con cable de poder (xyz) = N01 con cable de poder	○	5 VDC (D)
A black power adapter with a power cord.	As alternative: Medical CE approved power supply: Model: ME10A0503(xyz) (xyz) = B01 changeable direct plug-in (xyz) = Q01 with power cord (xyz) = N01 with power cord	●	5 VDC (1)

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