EARGRIPZ NO SLIPS

HEARING AID GRIPS
Double-sided adhesive pads for BTE and RIC hearing aids

FOR AUDIOLOGISTS AND DISPENSERS

Help patients keep their hearing aids on



For patients whose aids still slip after a good fitting. EarGripz are soft, double-sided hearing aid grips that help keep BTE and RIC aids on the ear through glasses, masks, sweaters and sweat.

No clips, no cords and nothing tied to the glasses

Made for the cases a fit change can't fix

Glasses wearers
Glasses go on and off all day, and the aid can go with them.

Mask wearers
Ear loops can flip the aid off when the mask comes off.

Soft or small ears
Soft cartilage, small ears, or ears that changed with age.

Memory care
A daily pad family or staff can put on, with no clip to lose.

Peel, stick, secure

Each pad holds in three places: on the hearing aid, behind the ear and against the side of the head.



1 Lift one pad

Peel back a corner of the top sheet. Lift one pad by its edges.



2 Stick it on

Press the center of the pad onto the side of the hearing aid.



3 Fold the wings

One wing to the back of the ear, one to the side of the head.

Helps the aid stay put

Clips and lanyards catch a hearing aid after it falls. EarGripz help hold it right where it sits, with nothing attached to clothing or glasses.

Best used alongside a good fitting, not in place of one.

PACK	DAILY USE	DESIGNED FOR	PAD SIZE	PACKAGING
60 pads, 2 sheets	1 pad per aid	BTE and RIC aids	20 mm (0.79 in)	Resealable zip pouch

Ask your distributor for EarGripz

Give patients a pack at the fitting, or stock it at the front desk.

Pehratek
855-EAR-8440 (855-327-8440)
info@pehratek.com · pehratek.com

Hearing aids not included. Use a new pad each day. Best hold on clean, dry skin.