



**Audiologist's Choice
AudioWipes**

ENVIRONMENTAL SURFACE WIPE

Contains LOW pH Synergized Quaternary Ammonium Chloride

Contains/CAS Number: Water/7732-18-5; Surfactant/68991-48-0; Citric Acid/77-92-9; Sodium Benzoate/532-32-1; _Alkyl dimethyl ethylbenzyl ammonium Chloride (C12-18)/68391-01-f; _Alkyl Dimethyl Benzyl Ammonium Chloride (C12-18)/68956-79-6

NFPA® Ratings: Health: 1; Flammability: 1; Reactivity: 0

DESIGNED FOR EFFECTIVE CLEANING:

Hearing aids, earmolds, headphones, ear inserts, specula, assistive listening devices, counter tops, telephones, instruments, sinks, toilets or any hard surface or object.

Safe to Use on All Contact Surfaces: Won't harm plastic or rubber

DIRECTIONS: Wipe desired surface with towel. Allow to air dry. No rinsing necessary.

DISPOSAL: Do not flush wipes

ITEM: #01301, 36 pop-up towelettes

MANUFACTURER: Oaktree Products, Inc
Chesterfield, MO 63005
800.347.1960

For additional information, see Safety Data Sheet