

Aurical® Otocam 300-III

Clarity in
every image



otometrics

a division of natus sensory

Video otoscope

A redesigned video otoscope with enhanced imaging and intuitive handling.



Aurical® Otocam 300-III



Greater image clarity

The new image system delivers clearer, wider visualization of the ear canal for more effective examination, counselling, and documentation.

Built for everyday durability

A reinforced cable and cable connection provide greater durability for busy clinical environments.

Refined ergonomics

New tactile buttons improve comfort and control while helping reduce clinician fatigue during everyday examinations.

Built-in electrical isolation

Supports the latest medical electrical safety standards for safe, reliable operation without additional electrical isolation.

Refined through clinical experience

The latest Aurical® Otocam 300-III builds on the trusted Otocam platform with enhanced image clarity, refined ergonomics, greater durability, and built-in electrical isolation.

Designed for today's hearing care workflow

The Aurical® Otocam 300-III makes it easy to capture, view, and share high-resolution images of the ear canal throughout the patient journey. With enhanced image clarity and wider visualization of the ear canal, clinicians can document ear health, support patient counseling, and perform quick visual inspections during hearing aid fittings with greater ease.

Building on the trusted Otocam platform, the latest generation introduces refined ergonomics, a reinforced cable and cable connection, built-in electrical isolation, and clear images from the start with no warm-up fogging—delivering an even better everyday experience while preserving the familiar workflow, software integration, and compatibility with existing Otocam speculas.




For diagnostics and counseling

- Enhanced image clarity and wider visualization of the ear canal
- Refined ergonomics with new tactile buttons for greater comfort and control
- Reinforced cable and cable connection for everyday durability
- Built-in electrical isolation supporting the latest medical electrical safety standards
- Clear images from the start with no warm-up fogging
- Lightweight USB-powered, plug-and-play design
- Compatible with the latest version of Otosuite® and Aurical® Freefit
- Compatible with existing Otocam speculas
- No special drivers required



**Comfortable to hold.
Easy to control.**

Compatible with your existing speculas.




Integrated where you work

Aurical® Otocam 300-III integrates directly with Aurical® Freefit and Otosuite®, allowing clinicians to perform quick visual inspections, document findings and support probe tube placement without interrupting the workflow. It can also be accessed directly from Aurical Freefit for probe tube placement and quick visual inspections during PMM.

Part of the Aurical® ecosystem

Fully integrated with Otosuite®



The Aurical Otocam 300-III integrates seamlessly with Otosuite, bringing otoscopy into the same connected workflow used for hearing assessment, verification, counseling, and documentation. Capture, review, and share images without leaving your clinical workflow.



Scan the QR code
to learn more about
Aurical Otocam 300-III.

Confidence in every sense

Natus Sensory is a global leader in sensory care, supporting clinicians across hearing, balance, vision, and newborn care. Guided by a belief that meaningful moments begin with a sense, we create trusted solutions that enable confident clinical decisions and improve patient outcomes – from giving every baby the best possible start in life to helping millions of people reconnect with the world through hearing and balance care.

natussensory.com